




Antibody Array

#Multi-Screening #동시분석 #시간절약 #가격절약 #Sample 절약 #자동분석



최소 샘플, 최대 분석 Screening Tool

📄 동시 분석 최적화

동시분석에 최적화 된 Antibody pairs 사용
Analytes간의 간섭 없는 최상의 결과 제공

📉 최고의 경제적 선택

최소 50 μ l의 Sample 사용
단 한번의 실험으로 최대 119개 Analytes 확인

💻 전용 분석 Software

분석도 쉽고 빠르게, 클릭 몇 번으로 분석 완료


🚫 특별한 장비 불필요

Western blot 장비 (Chemi Doc, X-ray film)로 Detection


📄 월등한 Reference 수

👤 웅비메디텍 분석서비스


Sample만 맡기면 분석까지 한번에!
Analytes가 많아 번거롭다면 Reading Service만 선택적으로!




#하나의 Membrane에서 다수의 결과 확인 가능



C. Untreated vs IL-1α + OSM Treated



#실험도 빠르게, 결과분석은 더 빠르게 Quick Spot Tool



Quick Spot Tool 분석 화면

- ✓ Cat# 선택하고 위치만 지정해주면 끝!
- ✓ Spot 범위 지정, Analytes 이름 자동 입력!
- ✓ 상황에 맞는 3가지 구매 옵션 (60일, 1년, 영구소장)

#제품 LIST 확인하기



Human						
Related research	Description	# of Analyte	Cat #	Description	# of Analyte	Cat #
Cancer	Phospho-RTK	49	ARY001B	NFκB Pathway	45(41/4)	ARY029
	Phospho-Kinase	39(37/2)	ARY003C	Angiogenesis	55	ARY007
	XL Oncology	84	ARY026			
Inflammation	Cytokine	36	ARY005B	XL Cytokine	105	ARY022B
	Chemokine	31	ARY017	Phospho-Immuno R	59	ARY004B
Stem cell	Pluripotent SC	15	ARY010			
Diabete/Obesity	Adipokine	58	ARY024			
기타	Kidney Biomarker	38	ARY019	Apoptosis	35	ARY009
	Cell Stress	26	ARY018	Protease	35	ARY021B
	Protease Inhibitor	32	ARY023	Protease/Inhibitor	67(35/32)	ARY025
	Ubiquitin	49	ARY027	Non-hematopoietic sR	119(62/57)	ARY012

Mouse						
Related research	Description	# of Analyte	Cat #	Description	# of Analyte	Cat #
Cancer	Phospho-RTK	39	ARY014	Angiogenesis	53	ARY015
Inflammation	Cytokine	40	ARY006	XL Cytokine	111	ARY028
	Chemokine	25	ARY020	Apoptosis	21	ARY031
Diabete/Obesity	Adipokine	38	ARY013			

Rat						
Related research	Description	# of Analyte	Cat #	Description	# of Analyte	Cat #
Inflammation	Cytokine	29	ARY008	XL Cytokine	79	ARY030
Diabete/Obesity	Adipokine	30	ARY016			