

Bio-Techne 주요 Antibody Overview

연구 분야 별 필요한 **Antibody**를 찾아보세요.

Bio-Techne의 검증된 고품질 Antibody


- ▶ 신뢰할 수 있는 품질로 연구 재현성 향상

Western Blot, Flow Cytometry, Immunostaining 등 폭넓은 application 지원

- ▶ 하나의 antibody로 다양한 연구에 적용 가능




Western Blot




Flow Cytometry



IHC/ICC



Neutralization



Matched Pairs

Cancer

Target	Cat#	Applications
*E-Cadherin	MAB18381	WB, Flow, IHC, CyTOF, ICC
*†Ras	NBP3-15909	WB, ELISA, ICC/IF
*†Ras - G12D Mutant	NBP3-13647	WB, IHC, ICC/IF
*Ras - G12V Mutant	NBP3-13668	WB, IHC
*EGFR	NBP3-32255	WB, Flow, IHC, ICC/IF
*†CDK4	NB120-6315	WB, IHC, ICC/IF, IP, MA
*†PIK3CA	NBP3-26384	WB, ICC/IF, ELISA
*†Ki-67/MKI67	NBP2-54791	WB, Flow, IHC, ICC/IF, KO
*†p53	NB200-103	WB, Flow, IHC, ICC/IF, ELISA, IP, CyTOF
*ErbB2/Her2	MAB1129	WB, Flow, IHC, ICC, CyTOF
*†VEGF	NB100-664	WB, Flow, IHC, ICC/IF, ELISA

* Human, † Mouse

Immunology

Target	Cat#	Applications
*NKp46/NCR1	MAB1850	WB, Flow, IHC, ICC, ISH, AgAct, CyTOF
*†CD31/PECAM-1	NB600-562	WB, Flow, IHC, ICC/IF, ELISA, ISH
*†DPPIV/CD26	NBP2-78791	WB, Flow, IHC, ICC/IF, KO
*PD-L1/B7-H1	MAB1562	WB, Flow, Block, CyTOF, Neut
*PD-1	MAB10861	WB, Flow, Block, CyTOF
*CD27/TNFRSF7	NBP2-78837	WB, Flow, IHC, ICC/IF, IP
*CD39/ENTPD1	NBP2-50387	WB, Flow, ICC/IF, IP, CyTOF
*CD25/IL-2R	MAB1020	WB, Flow, Block, CyTOF
*CD163	NBP1-60158	WB, Flow, IHC, ICC/IF, ELISA, IP
*BAFF/BLYS/TNFSF13B	MAB1241	WB, Flow, CyTOF
*BAFFR/TNFRSF13C	NBP3-27127	WB, Flow, ELISA

* Human, † Mouse

Inflammation

Target	Cat#	Applications
*IL-2	MAB104421	WB, ICC, ELISA
*IL-6	MAB206	WB, IHC, Neut, Pair
*IL-17	MAB3171	WB, ICFlow, IHC, IP, CyTOF
*IL-22	NBP2-27339	WB, IHC
*†Oncostatin M(OSM)	NBP3-16686	WB, IHC, ICC/IF, ELISA
*TNFSF13/APRIL	NBP1-97587	WB, Flow, IHC, ICC/IF, IP
*Granzyme K	MAB10216	WB, Flow, IHC
*Erythropoietin(EPO)	NBP2-53392	WB, IHC, ELISA, Functional
*†TNF-alpha	NBP2-27223	WB, Flow, Functional
*†COX-2	NBP3-16220	WB, IHC, ELISA
*†NFkB p65	NB100-56712	WB, Flow, IHC, IP, CyTOF
*NFkB p105/p50	NBP2-22178	WB, Flow, IHC, ICC/IF, ELISA

* Human, † Mouse




더 많은 antibody를
원하신다면!

다양한 Conjugated Antibody 솔루션

손쉬운 선택, 폭넓은 Conjugation 옵션

Excitation Laser (nm)	Conjugate	Excitation Max (nm)	Emission Max (nm)	Emission Color
UV (350)	Alexa Fluor® 350	346	442	
UV (350)	DyLight™ 350	353	432	
Violet (405)	DyLight™ 405	400	420	
Violet (405)	Alexa Fluor® 405	401	421	
Violet (405)	mFluor® Violet 450	406	445	
Violet (405)	mFluor® Violet 500	410	501	
Violet (405)	mFluor® Violet 610	421	613	
Blue (488)	PerCP	482	678	
Blue (488)	PE/Cy7	489,566	782	
Blue (488)	Alexa Fluor® 488	490	525	
Blue (488)	DyLight™ 488	493	518	
Blue (488)	FITC	498	519	
Blue (488)	PE	498,565	578	
Blue (488)	PE/Atto594	498,566	632	
Blue (488)	PE/Cy5.5	498,866	700	
Y-G (561)	Alexa Fluor® 532	532	554	
Y-G (561)	Janelia Fluor® 549	549	571	
Y-G (561)	DyLight™ 550	562	576	
Y-G (561)	Alexa Fluor® 594	590	617	
Y-G (561)	DyLight™ 594	593	618	
Red (633,640)	Janelia Fluor® 646	646	664	
Red (633,640)	APC	650	660	
Red (633,640)	DyLight™ 650	652	672	
Red (633,640)	Alexa Fluor® 647	653	669	
Red (633,640)	DyLight™ 680	692	712	
Red (633,640)	Alexa Fluor® 700	696	720	
Red (633,640)	Alexa Fluor® 750	749	775	
Red (633,640)	APC/Cy7	652	790	
	Biotin			
	HRP			

Y-G: yellow-green



형광 조합이 어려우신가요?
Novus Spectra Viewer를
이용해 보세요!